

deep_learning_motion_mask_segmentation




Figure 1: image from web

This repository is used to demonstrate the method published in : <http://link-to-paper.com>

Authors : Ludmilla Penarrubia, Nicolas Pinon, Emmanuel Roux, Eduardo Enrique Davila Serrano, Jean-Christophe Richard, Maciej Orkisz and David Sarrut.

This repository has several usages :

- 1- Use our trained model on the data we provide, as a proof of concept.
- 2- Test our trained model on your data, to get the motion mask segmentations on your data.
- 3- Train our model on your data and test it on your data

Pre-requisites and installations

- Make sure you have python3 installed
- Clone this repository on your machine and go in it:

```
cd deep_learning_motion_mask_segmentation/
```
- Create a virtual environments

```
python3 -m venv motion_mask_seg
```
- Activate the virtual environment

```
source motion_mask_seg/bin/activate
```

- Update pip3 repository and install dependencies listed in the requirements.txt

```
pip3 install --upgrade pip
python3 -m pip install -r requirements.txt
```

- Install Gatetools for preprocessing (optional)

```
pip3 install gatetools
```

Case 1 : Use *our* trained model on *our* showcase data

Run `python3 trained_model_on_showcase_data.py`

Motion mask as .mhd and .raw files will be located in : `results_showcase/`

Case 2 : Use *our* trained model on *your* data (work in progress)

Put all your .nii or .mgh or .mhd in the directory `data/`

(optional) Run : `gatetools/bin/gt_affine_transform -i input_data.mhd -o output_data.mhd --newspacing "2.0" --force_resample --adaptative -p "-1000.0"`

Edit the file `infer_motion_masks.py` (l. 62-63) with the path to your data and its size:

```
python          params.input_img_path = "./data/PATH_OF_THE_IMAGE_TO_SEGMENT.mhd"
# EDIT THIS LINE      params.input_size = [256, 256, 256] #
EDIT THIS LINE
```

Run : `python3 infer_motion_masks.py`

Motion mask as .nii files and figures will be located in : `results/inference_on_your_data/`

We suggest skipping the preprocessing step only if your data is sampled as isotropic 2mm³

(Advanced) Case 3 : Train and test our model on *your* data

For this use-case, we recommend plugging in your code the model located in `model.py`, we do not provide the data management part of the code, as it is really specific to each user.

Acknowledgments

Thanks to the authors of this repository : <https://github.com/milesial/Pytorch-UNet> for providing an efficient implementation of U-net.

This work was performed within the framework of the LABEX PRIMES (ANR-11-LABX-0063) of Université de Lyon, within the program “Investissements d’Avenir”(ANR-11-IDEX-0007) operated by the French National Research Agency (ANR).

Thanks to Olivier Bernard for getting us started with the project by providing examples of his codes.